

Production Order: 500000294239



Production Order Document
Production Order Qty: 500
PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Material Type:	ZFRT	Description: Edwards Flex Shaft Commander 155885	Order Type: ZSTD
Production Version:	7987		Project Phase:
Plant / Business Unit:	1213 / AC5		

NC - 29096

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials																							
50	KITTING3 Kitting Devices  Kitting Devices	<p>Kitting Devices Perform Order Kitting, Load Minor Mandrels, Dry Extrusions, and Cut FEP Record Time Extrusions Enter Dryer (Initial/Time/Date): <u>Am 68 6:00 PM 05 Jan 24</u> Record Time Extrusions First Exit Dryer (Initial/Time/Date): <u>Am 68 9:30 AM 06 Jan 24</u> Record Dryer Shelf #: <u>N/A</u></p> <hr/> <table border="1"> <thead> <tr> <th>Component Number</th> <th>Req'd Rev Rev Used</th> <th>UOM</th> <th>Qty.</th> <th>Batch No.</th> <th>Actual Qty Used</th> </tr> </thead> <tbody> <tr> <td>MM0179-01</td> <td>D <u>D</u></td> <td>PC</td> <td>500</td> <td><u>0000276172</u></td> <td><u>500</u></td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td><u>N/A</u></td> <td><u>N/A</u></td> </tr> <tr> <td>MM1536-01</td> <td>B <u>B</u></td> <td>PC</td> <td>500</td> <td><u>0000281412</u></td> <td><u>500</u></td> </tr> </tbody> </table>	Component Number	Req'd Rev Rev Used	UOM	Qty.	Batch No.	Actual Qty Used	MM0179-01	D <u>D</u>	PC	500	<u>0000276172</u>	<u>500</u>					<u>N/A</u>	<u>N/A</u>	MM1536-01	B <u>B</u>	PC	500	<u>0000281412</u>	<u>500</u>	N/A	N/A	4 JAN 24 JAT
Component Number	Req'd Rev Rev Used	UOM	Qty.	Batch No.	Actual Qty Used																								
MM0179-01	D <u>D</u>	PC	500	<u>0000276172</u>	<u>500</u>																								
				<u>N/A</u>	<u>N/A</u>																								
MM1536-01	B <u>B</u>	PC	500	<u>0000281412</u>	<u>500</u>																								

Notes: DA 2564, 2484

N/A

N/A

Date Printed: 04.01.2024 / 18:24:30

Am 68 06 Jan 24


Page: 1 of 18



SA0155-01

CREGANNA
MEDICAL
is part of





Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details						Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	RM0158-01	E	<u>E</u>	PC	200	<u>81054</u> <u>N/A</u>	<u>N/A</u> <u>100</u>			
		1000-1153-01	A	<u>A</u>	PC	594	<u>867411</u> <u>86781</u> <u>86780</u>	<u>N/A</u> <u>200</u> <u>200</u> <u>200</u>			
		1000-2053-01	A	<u>A</u>	PC	500	<u>0000287543</u> <u>000027880</u>	<u>400</u> <u>100</u>			
		MM1537-02	A	<u>A</u>	PC	500	<u>0000276175</u> <u>N/A</u>	<u>500</u> <u>N/A</u>			N/A
		TL0167-02	E	<u>E</u>	PC	70	<u>N/A</u> <u>N/A</u>	<u>Bulk</u> <u>Bulk</u>	<u>N/A</u> <u>N/A</u>		
		TL0165-05	J	<u>J</u>	PC	5	<u>N/A</u> <u>N/A</u>	<u>Bulk</u> <u>Bulk</u>			
		TL0165-03	J	<u>J</u>	PC	5	<u>N/A</u> <u>N/A</u>	<u>Bulk</u> <u>Bulk</u>			

Notes:

N/A

N/A

N/A

Date Printed: 04.01.2024 / 18:24:30

Page: 2 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details						Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials	
N/A	N/A	141967-01	02	02	PC	500	85502	525				
		RM7349-02	C	c	PC	543	82673	N/A				
		RM7348-01	C	c	PC	500	82883	582				
		RM4001-01	B	B	PC	125	82824	600				
		RM0607-01	D	D	PC	56	82824	N/A				
		RM0498-01	C	c	PC	500	0000287519	200				
		RM0009-04	I	I	PC	1	82971	N/A				
		RM0009-04	I	I	PC	1	8742187971	450				
								65				
Notes:		N/A										
		N/A										
		N/A										

Date Printed: 04.01.2024 / 18:24:30

①CB58 12Jan24 Page: 3 of 18



SA0155-01

CREGANNA
MEDICAL
is part of





Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details					Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	MM1538-01	A	<u>A</u>	PC	500	<u>N/A</u>	<u>Bulk</u>		
							<u>0000271052</u>	<u>500</u>		
		MM1537-01	A	<u>A</u>	PC	1000	<u>N/A</u>	<u>N/A</u>		
							<u>0000281413</u>	<u>1000</u>		
		MM0177-01	C	<u>C</u>	PC	500	<u>N/A</u>	<u>N/A</u>		
							<u>0000278966</u>	<u>500</u>		
		MM0180-01	E	<u>E</u>	PC	500	<u>N/A</u>	<u>N/A</u>		
							<u>0000275691</u>	<u>500</u>	<u>N/A</u>	<u>N/A</u>
		MM0178-01	E	<u>E</u>	PC	500	<u>N/A</u>	<u>N/A</u>		
							<u>0000271050</u>	<u>500</u>		
		MM0176-01	D	<u>D</u>	PC	500	<u>N/A</u>	<u>N/A</u>		
							<u>0000281411</u>	<u>500</u>		
		MM0074-01	G	<u>G</u>	PC	500	<u>N/A</u>	<u>40</u>		
							<u>0000271036</u>	<u>518</u>		
							<u>0000286923</u>	<u>26</u>		

Notes:

N/A

N/A

N/A

Date Printed: 04.01.2024 / 18:24:30

Page: 4 of 18



SA0155-01

CREGANNA
MEDICAL
is part of

Production Order: 500000294239

Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	N/A	N/A	N/A	N/A	N/A
100	CATASY01 Catheter Assembly 1 	Line Clearance Perform Line Clearance and Heat Gun Setting	500	0	07Jan24 CB58	
	Line Clearance Confirmation Reqd(Milestone)					
150	CATASY01 Catheter Assembly 1 	Major and Minor Mandrel Assembly	500	0	01Jan24	SD34 YKHO <i>BLW</i>
	Major and Minor Mandrel Assembly					
Notes:						
N/A						
N/A						
N/A						

Date Printed: 04.01.2024 / 18:24:30

Page: 5 of 18



SA0155-01

CREGANNA MEDICAL
is part of

Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
200	CATASY01 Catheter Assembly 1 	Loading Braid Stock	500	0	07Jan24	MC17 SC10
	Loading Braid Stock					
	Confirmation Reqd(Milestone)					
250	CATASY01 Catheter Assembly 1 	Trim Braid Wire at Proximal End	500	0	07Jan24	PY67 AI65 PZ22 TRN RL47
Notes:						
N/A						
N/A						
N/A						

Date Printed: 04.01.2024 / 18:24:30

Page: 6 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Trim Braid Wire at Proximal End Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
300	CATASY01 Catheter Assembly 1 	Insert Cut Hypo Tube Insert Cut Hypo Tube Confirmation Reqd(Milestone)	500	0	07Jan24	ALL2 RL47 PL22 TRV BD64
350	CATASY01 Catheter Assembly 1	Load Tubing	500	0	07Jan24	BD64 CX65

Notes:

N/A

N/A

N/A

Date Printed: 04.01.2024 / 18:24:30

Page: 7 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Load Tubing Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
400	CATASY01 Catheter Assembly 1 Reflow Confirmation Reqd(Milestone)	Reflow	500	0	07Jan24	AL67 CD19 AK02 ① SN67
450	CATASY01 Catheter	FEP Removal	500	0	07Jan24	AM47 JC92 SD34
Notes:						
N/A						
N/A						
N/A						

Date Printed: 04.01.2024 / 18:24:30

① AK02 06 Jan 24

Page: 8 of 18



SA0155-01

CREGANNA
MEDICAL
is part of





Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	Assembly 1 					YK40
N/A	FEP Removal	N/A	N/A	N/A	N/A	
	Confirmation Reqd(Milestone)					
500	CATASY01 Catheter Assembly 1 	In-process Inspection and Rework Material Consumed: Part #: 1000-1153-01 Batch #: 86741 Qty: N/A Part #: N/A Batch #: N/A Qty: N/A	DF75	FM-111 DH-111 OF-144-11 EW-144-14 15 25	PT09 AR02 BI60 KK54 YK95 08 Jan 24	
N/A	In-process Inspection and Rework					
	Confirmation Reqd(Milestone)					
N/A	N/A	N/A	N/A	N/A	N/A	N/A
Notes:						
N/A						
N/A						
N/A						

Date Printed: 04.01.2024 / 18:24:30

Page: 9 of 18



SA0155-01

CREGANNA MEDICAL
is part of

Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
550	CATASY01 Catheter Assembly 1 	Remove Heat Shrink & Mandrel Remove Heat Shrink & Mandrel Confirmation Reqd(Milestone)	468	MAH-HCI IDB-1 ⑨	06 Jun 24	PT09 MH10 CB19 TDN SGG MH10 YK95
600	CATASY01 Catheter Assembly 1 Distal Tip Assembly Confirmation	Distal Tip Assembly	468	○	06 Jun 24	PT09

Notes:

N/A

N/A

N/A

Date Printed: 04.01.2024 / 18:24:30

Page: 10 of 18



SA0155-01

CREGANNA MEDICAL
is part of



Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
650	CATASY01 Catheter Assembly 1 	Loading Heat Shrink	468	0	08Jan24	PT09 LH45
	Loading Heat Shrink					
	Confirmation Reqd(Milestone)					
700	CATASY01 Catheter Assembly 1 	Tipping Record Tipping Oven Information: TMI: <u>0386</u> Cal Due: <u>31may24</u> TMI: <u>0936A</u> Cal Due: <u>31may24</u> TMI: <u>0521</u> Cal Due: <u>31may24</u> TMI: <u>2083C</u> Cal Due: <u>31may24</u> Tipping	468	0	08Jan24	YH40 Hv36
Notes:						
N/A						
N/A						
N/A						

Date Printed: 04.01.2024 / 18:24:30

Page: 11 of 18



SA0155-01

CREGANNA
MEDICAL
is part of





Material: SA0155-01 Rev F


Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
750	CATASY01 Catheter Assembly 1 	Tip Inspection/ Flash Removal Material Consumed: Part # RM0607-01 Batch #: 78316 Qty: N/A Part #: RM4001-01 Batch #: 82824 Qty: N/A Part #: N/A Batch #: N/A Qty: N/A Part #: N/A Batch #: N/A Qty: N/A Part #: N/A Batch #: N/A Qty: N/A	468	0	08Jan24	PT09 BI60 Hv36 DY29
800	CATASY01 Catheter Assembly 1 	Major Mandrel Removal	459	ADC-#44111 ⑦①⑨	08Jan24	KL45 TRN SNL67 YK40 AI65 SS44
Notes:						
N/A						
N/A						
N/A						

Date Printed: 04.01.2024 / 18:24:30

Page: 12 of 18



SA0155-01

CREGANNA MEDICAL
is part of

(1) KL45 08Jan24



Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Major Mandrel Removal Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
850	CATASY01 Catheter Assembly 1  Cut to Length Confirmation Reqd(Milestone)	Cut to Length Record DIM05 gage result for the first 5 parts at the start of operation: 1. Pass 2. Pass 3. Pass 4. Pass 5. Pass	459	O	08 Jan 2024	RL45 AL42 TRN SS52
900	QUALITY1 Quality Inspection & Review	Quality Inspection and Review Perform Quality Inspection per QIP Document #3107610 Record Data in SAP ROS	N/A	N/A	N/A	N/A
Notes:						
N/A						
N/A						
N/A						

Date Printed: 04.01.2024 / 18:24:30

Page: 13 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Quality Inspection & Review Confirmation Reqd(Milestone)	<p>Re-Inspect after re-work.</p> <p>Required Inspection Visual/OD Inspection Record Inspection Data in SAP ROS Record Laser Micrometer Information: TMI: 0700-01 Cal Due: 31 May 24 TMI: N/A Cal Due: N/A TMI: N/A Cal Due: N/A Material Consumed: Part #: RM0607-01 Batch #: 78316 Qty: N/A Part #: RM4001-01 Batch #: 82824 Qty: N/A Part #: RM0158-01 Batch #: 81054 Qty: N/A Part #: 1000-153-01 Batch #: 86781 Qty: N/A Part #: N/A Batch #: N/A Qty: N/A</p>	448	SCR-1 EW-1 DIS-HH1 Fn-11 DL-1 ⑪	W35 LS46 YK95 P446 24 85 ^{con}	
950	QUALITY1 Quality Inspection & Review	<p>Quality Inspection & Review Borescope Inspection Record Inspection Data in SAP ROS Record Tip Gage Information: TMI: N/A Cal Due: N/A Record Caliper Information:</p>	N/A	N/A	N/A	N/A

Notes:

N/A

N/A

N/A

Date Printed: 04.01.2024 / 18:24:30

Page: 14 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	 Quality Inspection & Review	TMI: N/A Cal Due: N/A Record DIM02 Go/No-Go Gage Information: TMI: 0691 Cal Due: 30Sep25 TMI: 0692 Cal Due: 30Sep25 Record DIM02 Inspection Results N = 54: Pass: 54 Fail: 0	409	#90S11 DIS-HH H(CSP) STR-HH HH11 DIS-HH(HH) HH11 (39)	08Jun24	1267 Q11b P122 P146
1000	 Quality Inspection & Review	Quality Inspection & Review Leak Test Record Inspection Data in SAP ROS Record Leak Tester Information: TMI: 1056 Cal Due: 31May24 Record Length Gage Information: TMI: 0889D Cal Due: 30Sep24 Record Calibrated Ruler Information: TMI: 0629 Cal Due: 30Sep24	384	LT-HH HH H SCR-HH (TT) DEL-1// (TT) FM-11 (TT) VD-1 (25)	08Jun24	CB58

Notes:

N/A

N/A

N/A



SA0155-01



Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	N/A	N/A	N/A	N/A	N/A
1050	QUALITY1 Quality Inspection & Review 	Required Inspection Visual Final Inspection Perform Quality Inspection per QIP Document #3107610 Record Data in SAP ROS Confirmation Reqd(Milestone)	① 370 340	<ul style="list-style-type: none"> • FB - HII • EW - I • MEX - III • SKV - II • SCR - VHHII • VD - IIII • FM - I • DL - IIII • Del - II ① 144 • #50S - I • DIS - I • DT - I • STN - I • FL - I • GMR GVI - I <p>08 Jan 24</p>	SV43 XN26 KWS	
1100	CATASY01 Catheter Assembly 1 Line Closure	Line Closure Perform Line Closure Settle materials issued to production order (Initials/Date): GS05 08 Jan 24	N/A	N/A	08 Jan 24 GS05	
Notes:						
N/A						
N/A						
N/A						

Date Printed: 04.01.2024 / 18:24:30

Page: 16 of 18

① SV43 09 Jan 24
② SV43 09 Jan 24



SA0155-01

CREGANNA MEDICAL
is part of

Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
1150	PACKINT1 Packing assembly Package Confirmation Reqd(Milestone)	Package Package, Label, and Ship Finished Parts	346	O 10 Jan 24	AP10	

Notes:

N/A AP10 10 Jan 24

Date Printed: 04.01.2024 / 18:24:30

Page: 17 of 18



SA0155-01

CREGANNA MEDICAL
is part of



Production Order: 500000294239



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Batch Number: 0000294239

By: AP10

Date: 10 Jan 24

Reviewed By:

RB29

Date:

12 JAN 24

Notes:

N/A AP10 10 Jan 24/

Date Printed: 04.01.2024 / 18:24:30

Page: 18 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Beta to 2024 2228 1/1/23



CREGANNA
MEDICAL
is part of
TE

Extend to 2023 2228 Notes:

Beta to 2023 2228 1/1/23

Requestor Name: Udhesh Kapadnis

CONTROLLED COPY DEVIATION AUTHORIZATION NUMBER: 2484

* See attached email extension to 24 SEP 2023

TS2
24 AUG 23
23 OCT 2023 2228 1/1/23

J228

24 SEP 23

J228

25 APR 23

J228

24 APR 23

J228

J228

J228

J228

J228

J228

Document Number Affected	Revision
3107610	L

Deviation From:
QIP3107610, Section 8.0 Inspection Requirements (Supplemental Visual Inspection) OP 1050:
Current QIP3107610 does not state to inspect for the correct extrusion configuration.

Deviation To:
This DA allows addition inspection for correct assembly of extrusion material MM0179-01 and MM1536-01 during Performing QIP3107610, Section 8.0 Inspection Requirements (Supplemental Visual Inspection) OP 1050. See instructions attached to this DA.

Pre-Release Assembly Diagram

Justification: Recently it has been found that operators are incorrectly assembling MM0179-01 and MM1536-01. The event documents in NC-26390, and NC-26426. Only few of experienced inspectors can detect finished unit that contains incorrect extrusion configuration, and inexperienced inspectors may not which potential non-conformance unit sent to customer. Interim correction action has been implemented at OP 250, 300, 350 to detect unit built with out of oriented extrusions. This DA is adding another layer of inspection at final QC inspection to avoid incorrect assembly defects.

Part Number Affected	Revision
SA0155-01	H

Start Date: 26 Jul 2023 **End Date:** 25 Aug 2023 **Lot Number:** N/A

Risk Assessment:

Is there any potential risk(s) that may occur as a result of the proposed deviation including the following:
Control Plans Yes No FMEA's Yes No Validations Yes No
Details (if any): N/A

If yes to any of the above, what controls are being put in place to mitigate the risk.

Corrective Action Required: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	If no, explain:
If no, explain:	

If no, explain: No corrective action is required for this event as there are no changes to the current process, consumption of material, or how the product is produced. This added inspection guidelines are to avoid incorrect extrusion assembly defects.

Training Required: Yes No **If no, explain:**

① UK55, 23JW 2023



DA	2484 2468
①	

Description/Objectives of Training:
DA- Inspection at final QC, Op#1050.

Procedure:

- 100% inspection at Op#1050 per the instructions below.
- Inspect 1 part at a time.
- Inspection is focused on the correct MM0179-01 and MM1536-01 assembly.
- Use the example MM0179-01 and MM1536-02 fixture for inspection. (See image 1)
① MM01536-01 type Correction TS12 10AUG23



Image- 1

Step 1:

- Visually locate the MM0180-01 (Vestamid) transition to MM0179-01 on the completed part approximately 9.75" from the distal end using magnification light 2.25X minimum.
- Align the fixture MM0179-01 extrusion proximal end to the Vestamid transition on completed part. (See image 2)



Image- 2

- Visually verify the MM0179-01 distal end of the fixture is approximately at the same location on the completed part. (See image 3)




Image- 3

- Scrap the part if the transition is not approximately aligned. Save the scrapped parts for Engineer review.
- If the part transition is aligned, move to Step 2.

CONTROLLED COPY

Step 2:

- Visually verify the MM1536-01 distal end of the fixture is approximately at the same location on the completed part. (See image 4)




Image- 4

- Scrap the part if the transition is not approximately aligned. Save the scrapped parts for Engineer review.
- If the part transition is aligned, the part passes inspection.
- Use Image 5 as a guide for GOOD and BAD extrusion transition alignment.

1	MM0179-01 GOOD PART	MM1536-01
2	MM1536-01	MM0179-01 MM0179-01 and MM1536-01 Wrong Order - BAD PART
3	MM0179-01 Two MM0179-01 - BAD PART	MM0179-01
4	MM1536-01 Two MM1536-01 - BAD PART	MM1536-01

Image - 5

Title	Approval Name	Approval Signature	Date
Mgr. Quality Engineering	Hai Nguyen		25 Jul 2023
Mgr. Manufacturing Engineering	Jake Stanislawski		25 Jul 2023
Mgr. Operations	Matthew Benson		25 Jul 2023

FM0002.RevF

Deviation Authorization

CONTROLLED COPY

CONTROLLED COPY

DEVIATION AUTHORIZATION FORM


Requestor Name: Krishna Selvaraj

Document Number Affected	Revision
Doc #3005206 (MPI0238)	BP
Deviation From:	Deviation To:
Doc #3005206 (Flex Commander MPI0238): OPER850.11: Using a laser micrometer, check the DIM06 outer diameter. Position the laser indicator as close to the distal edge as possible. Start the measurement, then slowly move the part through the laser micrometer until reaching the lower edge of the shoulder.	Doc #3005206 (Flex Commander MPI0238): OPER850.11: Using a laser micrometer at OPER900 (TMI0700-01) , check the DIM06 outer diameter. Position the laser indicator as close to the distal edge as possible. Start the measurement, then slowly move the part through the laser micrometer until reaching the lower edge of the shoulder.

Justification:

TMI0602 lasermic which is currently used in SA0155-01 Flex commander product at OPER850 for Dim 6 inspection has mechanical failure and confirmed as not usable.

TM10700-01 lasermic is used at OPER900 for 100% inspection for Dim 1, Dim 6 and Dim 9. Since TM10700-01 is already qualified to inspect Dim 6 per ES0647: Laser micrometer equivalency test, there is no additional risk in using TMI0700-01 for OPER850 Dim 6 inspection till TMI0602 issue is resolved.

Part Number Affected	Revision	
SA0155-01	H	
Start Date:	End Date:	Lot Number:
16 Nov 23	15 DEC 23	N/A
Risk Assessment: Is there any potential risk(s) that may occur as a result of the proposed deviation including the following: Control Plans <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No FMEA's <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Validations <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Details (if any): N/A		
If yes to any of the above, what controls are being put in place to mitigate the risk – N/A		
Corrective Action Required: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No		
If no, explain: This is a temporary change to use TMI0700-01. DA will be removed once the lasermic TMI0602 issues are resolved and accepted for usage.		
Training Required: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No If no, explain: N/A		
Title	Approval Name	Approval Signature
Engineering Manager	Jake Stanislowski	
Quality Manager	Jay Zabel	
Operations Manager	Matthew Benson	



Document No: 5105589
FM15104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 500000294239

OP 400

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
Tm10942	44	10:20am	430	cm99	06Jan24	10:32am	415	cm99	06Jan24	16
Tm10942	44	10:47am	430	cm99	06Jan24	10:57am	415	cm99	06Jan24	16
Tm10942	44	11:16am	430	cm99	06Jan24	11:28am	415	cm99	06Jan24	16
Tm10942	44	11:42am	430	cm99	06Jan24	11:54am	415	cm99	06Jan24	16
Tm10942	44	12:11pm	430	cm99	06Jan24	12:12pm	415	cm99	06Jan24	16
Tm10942	44	1:30 PM	430	RL47	06Jan24	1:42 PM	415	RL47	06Jan24	16
Tm10942	44	2:33PM	430	YKH0	06Jan24	2:45PM	415	YKH0	06Jan24	16
Tm10942	44	3:02 PM	430	RL47	06Jan24	3:14 PM	415	RL47	06Jan24	16
Tm10942	44	4:46PM	430	SD34	06Jan24	4:58pm	415	SD34	06Jan24	16
Tm10942	44	5:37PM	430	RL47	06Jan24	5:49PM	415	RL47	06Jan24	16
Tm10942	44	6:27PM	430	CD19	06Jan24	6:39PM	415	CD19	06Jan24	16
Tm10942	44	7:06am	430	cm99	07Jan24	7:18am	415	cm99	07Jan24	16



Document No: 5105589
FM15104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 50000294239

OP 400

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
Tm10942	44	8:16am	430	cmq99	07Jan24	8:28am	415	cmq99	07Jan24	16
Tm10942	44	9:41 AM	430	Rlb	07Jan24	9:53 AM	415	Rlb	07Jan24	16

cmq99 07Jan24
NIA cmq9



Document No: 5105589
FM5104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 500000294239

OP 400

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
TM10745	44	10:32am	430	cm99	06Jan24	10:44am	415	cm99	06Jan24	16
TM10745	44	11:01am	430	cm99	06Jan24	11:13am	415	cm99	06Jan24	16
TM10745	44	11:34am	430	cm99	06Jan24	11:46am	415	cm99	06Jan24	16
TM10745	44	11:58am	430	cm99	06Jan24	12:10pm	415	cm99	06Jan24	16
TM10745	44	12:23pm	480	SN67	06Jan24	12:38pm	415	SN67	06Jan24	16
TM10745	44	1:22pm	430	RL47	06Jan24	1:34pm	415	RL47	06Jan24	16
TM10745	44	1:46pm	430	RL47	06Jan24	1:58pm	415	RL47	06Jan24	16
TM10745	44	2:07pm	430	SD34	06Jan24	2:18pm	415	SD34	06Jan24	16
TM10745	44	2:47pm	438	SN67	06Jan24	2:59pm	415	SN67	06Jan24	16
TM10745	44	3:23pm	430	SN67	06Jan24	3:35pm	415	SN67	06Jan24	16
JM10745	44	4:24pm	430	RL47	06Jan24	4:36pm	415	RL47	06Jan24	16
TM10745	44	5:02pm	430	SD34	06Jan24	5:14pm	415	SD34	06Jan24	16



Document No: 5105589
FM5104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# SO0000294239

OP 400

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
TM10745	44	5:22pm	430	RL47	^{06 Jan 24} 06 Jan 23 ②	5:34pm	415	RL47	^{06 Jan 24} 06 Jan 23 ②	16
TM10745	44	5:47pm	430	AL42	^{06 JAN 23} 06 Jan 24 ①	5:59pm	415	AL42	^{06 JAN 23} 06 Jan 24 ②	16
TM10745	44	6:18pm	430	RL47	^{06 Jan 24} 06 Jan 23 ②	6:30 pm	415	RL47	^{06 Jan 24} 06 Jan 23 ②	16
TM10745	44	7:36AM	430	SN67	^{07 JAN 24} 07 Jan 23 ①	7:48AM	415	SN67	^{07 JAN 24} 07 Jan 23 ②	16
TM10745	44	8:37AM	430	SN67	07 Jan 24	8:49AM	415	SN67	07 Jan 24	16
TM10745	44	9:45AM	430	Pv16	07 Jan 24	9:57 AM	415	RJ16	07 Jan 24	4

① SN67 07 Jan 24

② PV 08 Jan 24 correction for ~~RL47, AL42, RL47~~ ②



Document No: 5106073
Rev: E
Document Type: Manufacturing Form
Title: SA0155-01 Visual Rework Form

PO #: 500000294239

OP #: 500 Shift #: 3

Total Parts Reworked:		205	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles		25
EH	Exposed Hypotube		40
EW	Exposed Wire		85
MP	Micropores	N/A	N/A
SCR	Scratch		18
SKV	Skive Marks		15
VD	Voids		22
N/A	N/A	N/A	N/A
Inspected By (Sign and Date):		LS46, AR02	07 Jan 24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only.

Data Uploaded for Engineering Review (Check):



Document No: 5106073
Rev: E
Document Type: Manufacturing Form
Title: SA0155-01 Visual Rework Form

PO #: 500000 294239 OP #: 500 Shift #: 3

Total Parts Reworked:		37	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	N/A	N/A
EH	Exposed Hypotube		4
EW	Exposed Wire	11	27
MP	Micropores	N/A	N/A
SCR	Scratch	11	2
SKV	Skive Marks	1	1
VD	Voids	11	3
N/A	N/A	N/A	N/A
Inspected By (Sign and Date):		AR02	06 Jan 24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only.

Data Uploaded for Engineering Review (Check):



Document No: 6102646

Rev: A

Document Type: Manufacturing Form

Title: SA0155-01 Tipping Rework Form

PO #: 500000294239

OP #: 750 Shift #: 1st

Total Parts Reworked:		21	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
DIM07 OS / WO	DIM07 Oversized (Window Open)	HH1 HH1	11
DIM07 US / WC	DIM07 Undersized (Window Closed)	N/A	N/A
EH	Exposed Hypotube	III	3
		HH1 II	7
Inspected By (Sign and Date):		DY29	08 Jan 24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):



Document No: 6102646

Rev: A

Document Type: Manufacturing Form

Title: SA0155-01 Tipping Rework Form

PO #: 50000294239

OP #: 750 Shift #: 3

Total Parts Reworked:		170	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
DIM07 OS / WO	DIM07 Oversized (Window Open)	 	67
DIM07 US / WC	DIM07 Undersized (Window Closed)		18
EH	Exposed Hypotube		54
GD/AZ	Glue Damage / Air Bubbles		31
Inspected By (Sign and Date):		BIGO	08 Jan 24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):

500000294239

PRODUCTION ORDER# 500000294239 ① KL95 08Jan24
 (1)

OP 800

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
TM10409	N/A	7:12AM	190°F	YK40	07Jan24	8:22AM	190°F	YK40	07Jan24	31
TM12036	N/A	8:08AM	190°F	YK40	07Jan24	9:18AM	190°F	XK40	07Jan24	38
TM10409	N/A	9:51AM	190°F	YK40	07Jan24	11:51AM	190°F	YK40	07Jan24	36
TM12036	N/A	10:39AM	190°F	YK40	07Jan24	11:49AM	190°F	YK40	07Jan24	35
TM10409	N/A	11:36AM	190°F	YK40	07Jan24	12:46PM	190°F	YK40	07Jan24	32
TM12036	N/A	12:12pm	190°F	YK40	07Jan24	1:22PM	190°F	YK40	07Jan24	30
TM10409	N/A	1:30PM	190°F	AEG5	07Jan24	2:40PM	190°F	AEG5	07Jan24	21
TM12036	N/A	2:30PM	190°F	AEG5	07Jan24	3:40PM	190°F	AEG5	07Jan24	30
TM10409	N/A	5:50PM	190°F	CD19	07Jan24	6:06PM	190°F	CD19	07Jan24	53
TM10409	N/A	4:50am	190°F	SS44	08Jan24	6:00am	190°F	SS44	08Jan24	91
TM10409	N/A	6:05am	190°F	SS44	08Jan24	7:15am	190°F	SS44	08Jan24	62
N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

① YK40 07Jan24



Document No: 6102619

Rev: B

Document Type: Manufacturing Form

Title: SA0155-01 Dimensional/Visual Rework Form

PO #: 50000294239

OP #: 900 Shift #: 3

Total Parts Reworked:		121	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles		2
EH	Exposed Hypotube		6
EW	Exposed Wire		34
MP	Micropores	N/A	N/A
SCR	Scratch		29
SKV	Skive Marks		26
VD	Voids		8
DIM01 US	DIM01 OD Undersized	N/A	N/A
DIM06 US	DIM06 OD Undersized		12
DIM06 OS	DIM06 OD Oversized		4
DIM09 US	DIM09 OD Undersized	N/A	N/A
Inspected By (Sign and Date):		LS46, KX54, VJ35	07JAN24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):



Document No: 6102619

Rev: B

Document Type: Manufacturing Form

Title: SA0155-01 Dimensional/Visual Rework Form

PO #: 500000294239

OP #: 900 Shift #: 1st

Total Parts Reworked:		136	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	n/a	n/a
EH	Exposed Hypotube		11
EW	Exposed Wire		50
MP	Micropores	n/a	n/a
SCR	Scratch		90
SKV	Skive Marks	n/a	n/a
VD	Voids		21
DIM01 US	DIM01 OD Undersized	n/a	n/a
DIM06 US	DIM06 OD Undersized	n/a	n/a
DIM06 OS	DIM06 OD Oversized	n/a	n/a
DIM09 US	DIM09 OD Undersized	n/a	n/a
Inspected By (Sign and Date):		P1460429, KT47 K155 08 Jan 24	

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):

Maximum Force Reached During Tensile Test (10 samples accepted from final inspection for each lot shall be selected and tensile tested)																
Sample # -->	1	2	3	4	5	6	7	8	9	10	Avg	St Dev	K	Calculated Lower bound	Min Spec	Pass / Fail
Seg A	25.15	25.73	24.98	24.01	26.7	27.2	25.75	24.74	24.67	25.35	25.428	0.960067	4.378	21.2248261	8.542	PASS
Seg B	59.83	56.67	63.62	53.84	62.86	56.16	57.83	61.65	62.66	56.94	59.206	3.377363	3.981	45.7607171	8.542	PASS
Seg C	83.9	79.61	76.19	78.53	80.09	83.94	83.44	78.56	80.47	77.25	80.198	2.766252	2.911	72.1454399	8.542	PASS

All Force Values are recorded in Pound-Force and Distance is in Inches
Specification for lower bound is 38N was converted to 8.542Lbf
First Peak Force was collected during test and has been included in the raw data file (this information will not be captured / summarized in the DA due to it is not required to used for DA acceptance.

EDW Commander Flex - Bend and Tensile Strength Testing

LOT #: 500000294239

Date: 09 JAN 24

Inspector Name: Javier Olivares

Equipment ID: TMI0311B

Cal Due Date: 27 OCT 24

J001 09 jan 24

TE MEDICAL REWORK/REINSPECTION WORK ORDER FORM

Customer: Edwards Lifesciences	PO No.: N/A
Description: Flex Commander	TE Medical Lot No.:294239
Customer P/N: 155885	Returned Quantity: 384
TE Medical P/N: SA0155-01	Other: N/A
RMA #:	NC or CAPA #: NC29096
Will the rework have adverse effects on the product? <input checked="" type="checkbox"/> No <input type="checkbox"/> Yes (Provide justification of rework and how it was concluded).	
<p>There will be no product rework allowed for this NC. There will be no adverse effect on the products as quality inspectors are certified to perform 100% visual inspection per existing Quality Inspection Plan (QIP3107610)</p>	

Special Instructions:

- QE/ Technician: to complete training per special instructions before sorting activity. Attach training record to NC before closing.
- Inspectors: to perform the line clearance per MPI0230 prior to re-inspection of each lot.
- Inspectors: to verify lot quantity prior to performing re-inspection of the affected lots. Contact Quality if count discrepancy found.
- Inspectors: to perform 100% dimensional inspection of Dim5 using calibrated Flex Shaft Tip Go/No Go gauge TM150713 and Vernier Caliper for DIM7. Defective units to be identified with defect type.
- DA Inspector: to perform 100% visual inspection using visual guideline doc #3107585, before moving good units to Packaging. Contact Quality to remove hold in SAP.
- QC Operator: Confirm lots in SAP to reflect final accept/reject quantities from 200% inspection activities.
- Accepted units can be sent to Packaging and onward to Shipping.

Note: The scrap unit(s) will be documented within this Form and may be used for Engineering study (Quality Engineer is to collect scrapped parts)

Process	Qty	Initials	Equipment Number and Cal Due Date if applicable	Comments
Dimension inspection	384	P146	① 44A TM150713	N/A
Reworked (if necessary)	N/A	N/A		N/A
AQL Sampling	N/A	N/A		N/A
Total accepted	384	P146		N/A
Total failed / scrapped	0	P146		N/A
Other	N/A	N/A		N/A

①P146 09 Jan 24

TE MEDICAL

REWORK/REINSPECTION WORK ORDER FORM

Note: Fill in the table below if components were replaced during the rework process. This form will be filed along with the original LHR.

Component(s) Used			Component(s) removed/scrapped		
P/N	Rev.	Lot No.	Qty.	P/N	Rev.
NA P146	01 Jan 24				

Process	Time	Hour(s)	Initials	Comments
Re-work hours (\pm .25 hrs)	.5	PM	N/A	N/A
Final Rework/Re-Inspection/Re-Test hours (\pm .25 hrs)	N/A	N/A	N/A	N/A

Re-work/Re-inspection/Re-Test Performed By: Py46 Date: 09 Jan 24
 Final Rework/Re-Inspection/Re-Test Reviewed By: N/A Date: N/A

Disposition of Reject Material (unit that does not meet the acceptance criteria after rework/inspection):

Scrap to Bin

Scrap Disposition by (Name, Date & Sign): Alfonso Bahin 09 Jan 2024 Q. Bahin

Scrap to R&D ENGINEERING

AB25 09 Jan 2024