

Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Material Type:	ZFRT	Description: Edwards Flex Shaft Commander 155885	Order Type: ZSTD
Production Version:	7988		Project Phase:
Plant / Business Unit:	1213 / AC5		

Opr No.	Planned WorkCenter Description	Operation Details						Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
		Component Number	Req'd Rev Rev Used	UOM	Qty.	Batch No.	Actual Qty Used				
50	KITTING3 Kitting Devices 	Kitting Devices Perform Order Kitting, Load Minor Mandrels, Dry Extrusions, and Cut FEP Record Time Extrusions Enter Dryer (Initial/Time/Date): <u>GS85 7:00 AM</u> 24 Jan 24 Record Time Extrusions First Exit Dryer (Initial/Time/Date): <u>XC31 11:30pm</u> 25 JAN 24 Record Dryer Shelf #: <u>N/A</u>						N/A	N/A	25 JAN 24	SHJ
		Component Number	Req'd Rev Rev Used	UOM	Qty.	Batch No.	Actual Qty Used				
		1000-2053-01	A <u>A</u>	PC	500	<u>0000279380</u>	<u>500</u>				
						<u>N/A</u>	<u>N/A</u>				
		MM1537-02	A <u>A</u>	PC	500	<u>0000288401</u>	<u>500</u>				

Notes: DA 2484, 2564

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 1 of 18



SA0155-01

CREGANNA MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500
PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details					Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	RM0158-01	E	<u>E</u>	PC	200	<u>N/A</u>	<u>N/A</u>			
						<u>58497</u>	<u>200</u>			
	TL0167-02	E	<u>E</u>	PC	70	<u>N/A</u>	<u>N/A</u>			
						<u>N/A</u>	<u>Bulk</u>			
	TL0165-05	J	<u>J</u>	PC	5	<u>N/A</u>	<u>Bulk</u>			
						<u>N/A</u>	<u>Bulk</u>			
	TL0165-03	J	<u>J</u>	PC	5	<u>N/A</u>	<u>Bulk</u>	<u>N/A</u>	<u>N/A</u>	<u>29Jan24</u>
						<u>N/A</u>	<u>Bulk</u>			<u>KP02</u>
N/A	141967-01	02	<u>02</u>	PC	500	<u>85793</u>	<u>525</u>			
	RM7349-02	C	<u>c</u>	PC	543	<u>82834</u>	<u>N/A</u>			
						<u>82836</u>	<u>200</u>			
						<u>82856</u>	<u>300</u>			
							<u>100</u>			
	RM7348-01	C	<u>c</u>	PC	500	<u>85428</u>	<u>500</u>			
						<u>N/A</u>	<u>N/A</u>			

Notes:

N/A

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 2 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details						Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	RM4001-01	B	<u>B</u>	PC	125	<u>82469</u>	<u>200</u>			
		RM0607-01	D	<u>D</u>	PC	56	<u>74663</u>	<u>N/A</u>	<u>N/A</u>		
		RM0498-01	C	<u>C</u>	PC	500	<u>0000275492</u>	<u>500</u>	<u>N/A</u>		
		RM0362-01	E	<u>E</u>	PC	594	<u>80236</u>	<u>600</u>	<u>N/A</u>		
		RM0009-04	I	<u>I</u>	PC	1	<u>79168</u>	<u>Bulk</u>	<u>N/A</u>	<u>N/A</u>	<u>N/A</u>
		RM0009-04	I	<u>I</u>	PC	1	<u>79168</u>	<u>Bulk</u>	<u>N/A</u>		
		MM1538-01	A	<u>A</u>	PC	500	<u>0000278970</u>	<u>500</u>	<u>N/A</u>		
		MM1537-01	A	<u>A</u>	PC	1000	<u>0000284209</u>	<u>1000</u>			

Notes:

N/A

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 3 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details						Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	MM1536-01	B	<u>B</u>	PC	500	<u>0000281412</u>	N/A	N/A		
		MM0180-01	E	<u>E</u>	PC	500	<u>0000287541</u>	N/A	N/A		
		MM0179-01	D	<u>D</u>	PC	500	<u>0000276172</u>	N/A	N/A		
		MM0178-01	E	<u>E</u>	PC	500	<u>0000276174</u>	N/A	N/A	N/A	N/A
		MM0177-01	C	<u>C</u>	PC	500	<u>0000284208</u>	N/A	N/A		
		MM0176-01	D	<u>D</u>	PC	500	<u>0000281411</u>	N/A	N/A		
		MM0074-01	G	<u>G</u>	PC	500	<u>0000293933</u>	N/A	N/A		

Notes:

N/A

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 4 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	N/A	N/A	N/A	N/A	N/A
100	CATASY01 Catheter Assembly 1 	Line Clearance Perform Line Clearance and Heat Gun Setting	500	0	26Jan24	V078 CB58
	Line Clearance Confirmation Reqd(Milestone)					
150	CATASY01 Catheter Assembly 1 	Major and Minor Mandrel Assembly	500	0	26Jan24	CL30 JY90 SD34 YK40
	Major and Minor Mandrel Assembly					
Notes:						
N/A						
N/A						
N/A						

Date Printed: 23.01.2024 / 19:46:14

Page: 5 of 18



SA0155-01

CREGANNA MEDICAL
is part of



Production Order: 500000297999Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
200	CATASY01 Catheter Assembly 1  Loading Braid Stock Confirmation Reqd(Milestone)	Loading Braid Stock	500	0	26Jan24	CLOS ⁵ Y04 CD19 MV50
250	CATASY01 Catheter Assembly 1  Trim Braid Wire at Proximal End		500	0	26Jan24	ePSZ DX35 PY67 AL34

Notes:

N/A

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 6 of 18



SA0155-01

CREGANNA MEDICAL
is part of

Production Order: 500000297999Production Order Document
Production Order Qty: 500PC
Sheet: 1 of 1**Material: SA0155-01 Rev F**

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	Trim Braid Wire at Proximal End N/A Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
300	CATASY01 Catheter Assembly 1 	Insert Cut Hypo Tube	500	0	26Jan24	SH23 G522 AIBS SCIU
	Insert Cut Hypo Tube Confirmation Reqd(Milestone)					
350	CATASY01 Catheter Assembly 1	Load Tubing	500	0	26Jan24	ST96 BD64 RL47
Notes:						
N/A						
N/A						
N/A						

Date Printed: 23.01.2024 / 19:46:14

Page: 7 of 18



SA0155-01

CREGANNA MEDICAL
is part of

Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	 Load Tubing Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
400	CATASY01 Catheter Assembly 1 Reflow Confirmation Reqd(Milestone)		500	0	26Jan24	SY47 V078 SNL67 AL67 AM47
450	CATASY01 Catheter	FEP Removal	500	0	26Jan24	SH85 AM47
Notes:						
N/A						
N/A						
N/A						

Date Printed: 23.01.2024 / 19:46:14

Page: 8 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	Assembly 1 					
N/A	FEP Removal	N/A	N/A	N/A	N/A	N/A
	Confirmation Reqd(Milestone)					
500	CATASY01 Catheter Assembly 1 	In-process Inspection and Rework Material Consumed: Part #: 1000-1153-01 Batch #: 87654 Qty: N/A Part #: N/A Batch #: N/A Qty: N/A	489	EW-1X111 OF-1 DS-1 DW-1 (11)	26Jan24 VL91 R66 TD45 LS46 LLG1 AR02	
N/A	In-process Inspection and Rework					
	Confirmation Reqd(Milestone)					
N/A	N/A	N/A	N/A	N/A	N/A	N/A
Notes:						
N/A						
N/A						

Date Printed: 23.01.2024 / 19:46:14

Page: 9 of 18



SA0155-01

CREGANNA MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
550	CATASY01 Catheter Assembly 1 	Remove Heat Shrink & Mandrel Remove Heat Shrink & Mandrel Confirmation Reqd(Milestone)	489	0	26Jan24	PP40 L445 TRN MH10 PH59
600	CATASY01 Catheter Assembly 1 Distal Tip Assembly Distal Tip Assembly Confirmation	Distal Tip Assembly	489	0	26Jan24	HT72 ML60 KT47 PH59

Notes:

N/A

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 10 of 18



SA0155-01

CREGANNA MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
650	CATASY01 Catheter Assembly 1 	Loading Heat Shrink	489	0	26Jan24	SV46 PT09 LH45
	Loading Heat Shrink					
	Confirmation Reqd(Milestone)					
700	CATASY01 Catheter Assembly 1 	Tipping Record Tipping Oven Information: TMI: <u>2083C</u> Cal Due: <u>31 May 24</u> TMI: <u>0936A</u> Cal Due: <u>31 May 24</u> TMI: <u>0386</u> Cal Due: <u>31 May 24</u> TMI: <u>0521</u> Cal Due: <u>31 May 24</u> Tipping	489	0	26Jan24	PH59
Notes:						
N/A						
N/A						
N/A						

Date Printed: 23.01.2024 / 19:46:14

Page: 11 of 18



SA0155-01

CREGANNA
MEDICAL
is part of





Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
750	CATASY01 Catheter Assembly 1 	Tip Inspection/ Flash Removal Material Consumed: Part #: RM4001-01 Batch #: 82469 Qty: N/A Part #: RMD607-01 Batch #: 74663 Qty: N/A Part #: N/A Batch #: N/A Qty: N/A Part #: N/A Batch #: N/A Qty: N/A Part #: N/A Batch #: N/A Qty: N/A	487	EH-11 ②	26Jan24	0429
800	CATASY01 Catheter Assembly 1 	Major Mandrel Removal	485	ACD-11 ②	26Jan24	PZ22 YKHO BDG4
Notes:						
N/A						
N/A						
N/A						

Date Printed: 23.01.2024 / 19:46:14

Page: 12 of 18



SA0155-01

CREGANNA
MEDICAL
is part of

Production Order: 500000297999Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	Major Mandrel Removal N/A Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
850	CATASY01 Catheter Assembly 1  Cut to Length Confirmation Reqd(Milestone)	Cut to Length Record DIM05 gage result for the first 5 parts at the start of operation: 1. <u>Pass</u> 2. <u>Pass</u> 3. <u>Pass</u> 4. <u>Pass</u> 5. <u>Pass</u>	485	0	27Jan24	DLOT
900	QUALITY1 Quality Inspection & Review	Quality Inspection and Review Perform Quality Inspection per QIP Document #3107610 Record Data in SAP ROS	N/A	N/A	N/A	N/A
Notes:						
N/A						
N/A						
N/A						

Date Printed: 23.01.2024 / 19:46:14

Page: 13 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Quality Inspection & Review Confirmation Reqd(Milestone)	<p>Re-Inspect after re-work.</p> <p>Required Inspection Visual/OD Inspection Record Inspection Data in SAP ROS Record Laser Micrometer Information: TMI: <u>0700-01</u> Cal Due: <u>31 MAY 2024</u> TMI: <u>N/A</u> Cal Due: <u>N/A</u> TMI: <u>N/A</u> Cal Due: <u>N/A</u> Material Consumed: Part #: <u>1000-1153-01</u> Batch #: <u>87654</u> Qty: <u>N/A</u> Part #: <u>RM0158-01</u> Batch #: <u>58497</u> Qty: <u>N/A</u> Part #: <u>RM4001-01</u> Batch #: <u>82469</u> Qty: <u>N/A</u> Part #: <u>RM0607-01</u> Batch #: <u>74663</u> Qty: <u>N/A</u> Part #: <u>N/A</u> Batch #: <u>N/A</u> Qty: <u>N/A</u></p>	461	(24)	27Jan24	P146 mc17 KX54 DL07
950	QUALITY1 Quality Inspection & Review	<p>Quality Inspection & Review Borescope Inspection Record Inspection Data in SAP ROS Record Tip Gage Information: TMI: <u>N/A</u> Cal Due: <u>N/A</u> Record Caliper Information:</p>	N/A	N/A	N/A	N/A

Notes:

N/A

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 14 of 18



SA0155-01

CREGANNA MEDICAL
is part of

Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	 Quality Inspection & Review Confirmation Reqd(Milestone)	TMI: <u>N/A</u> Cal Due: <u>N/A</u> Record DIM02 Go/No-Go Gage Information: TMI: <u>0691</u> Cal Due: <u>30Sep2025</u> TMI: <u>0692</u> Cal Due: <u>30Sep2025</u> Record DIM02 Inspection Results N = 54: Pass: <u>54</u> Fail: <u>0</u>	454	DIS-IH STR-II (7)	27Jan24	QW10
1000	 Quality Inspection & Review Quality Inspection & Review Confirmation Reqd(Milestone)	Quality Inspection & Review Leak Test Record Inspection Data in SAP ROS Record Leak Tester Information: TMI: <u>1056</u> Cal Due: <u>31MAY2024</u> Record Length Gage Information: TMI: <u>0889</u> Cal Due: <u>30Sep2024</u> Record Calibrated Ruler Information: TMI: <u>0629</u> Cal Due: <u>30Sep2024</u>	450	DAL-III LT-1 (4)	27Jan24	QW10 CB58

Notes:

N/A

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 15 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	N/A	N/A	N/A	N/A	N/A
1050	QUALITY1 Quality Inspection & Review  Quality Inspection & Review Confirmation Reqd(Milestone)	Required Inspection Visual Final Inspection Perform Quality Inspection per QIP Document #3107610 Record Data in SAP ROS	427	• Del-1111 • SCR-111 • DL-11 • DIS-11 • DNT-11 • EW-11 • BP-11 • SKV-11 • FL-1 • CRK-1 • Mex-1 • VD-1	29Jan24 (23)	SV43 XN26
1100	CATASY01 Catheter Assembly 1  Line Closure	Line Closure Perform Line Closure Settle materials issued to production order (Initials/Date): KPD2 29Jan24	N/A	N/A	29Jan24 KPD2	

Notes:

N/A

N/A

N/A

Date Printed: 23.01.2024 / 19:46:14

Page: 16 of 18



SA0155-01

CREGANNA MEDICAL
is part of



Production Order: 500000297999



Production Order Document

Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
1100	Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
1150	PACKINT1 Packing assembly Package Confirmation Reqd(Milestone)	Package Package, Label, and Ship Finished Parts	U27	0	29 Jan 24	APW

Notes:

N/A APW 29 Jan 24 /

Date Printed: 23.01.2024 / 19:46:14

Page: 17 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000297999



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Batch Number: 0000297999

By: AP10

Date: 29 Jan 24

Reviewed By:

RB29

Date:

30 JAN 24

Notes:

N/A AP10 29 Jan 24

Date Printed: 23.01.2024 / 19:46:14

Page: 18 of 18



SA0155-01

**CREGANNA
MEDICAL**
is part of



Batches to 2024 3228 U/A/23
Expiry to 19 Feb 2024 3228 V/A/23
CREGANNA MEDICAL is part of



EX-2023-3228-V-A/23

Extend to 22 Nov 2023 3228 P00223
Batches to 2023 3228 V/A/23

Requestor Name: Udhesh Kapadia

DEVIATION AUTHORIZATION NUMBER: 2484
** See attached email extension to 24 SEP 23*
TS12
24 AUG 23
23 OCT 2023 3228 U/A/23

Document Number Affected	Revision
3107610	L

Deviation From:
QIP3107610, Section 8.0 Inspection Requirements (Supplemental Visual Inspection) OP 1050:
Current QIP3107610 does not state to inspect for the correct extrusion configuration.

Deviation To:
This DA allows addition inspection for correct assembly of extrusion material MM0179-01 and MM1536-01 during performing QIP3107610, Section 8.0 Inspection Requirements (Supplemental Visual Inspection) OP 1050.
See instructions attached to this DA.

Pre-Relow Assembly Diagram:

Justification: Recently it has been found that operators are incorrectly assembling MM0179-01 and MM1536-01. The event documents in NC-26390, and NC-26426. Only few of experienced inspectors can detect finished unit that contains incorrect extrusion configuration, and inexperienced inspectors may not which potential non-conformance unit sent to customer. Interim correction action has been implemented at OP 250, 300, 350 to detect unit built with out of oriented extrusions. This DA is adding another layer of inspection at final QC inspection to avoid incorrect assembly defects.

Part Number Affected	Revision		
SA0155-01	H		
Start Date:	End Date:	Lot Number:	
26 Jul 2023	25 Aug 2023	N/A	

Risk Assessment: Is there any potential risk(s) that may occur as a result of the proposed deviation including the following: Control Plans <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No FMEA's <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Validations <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Details (if any): N/A	Corrective Action Required: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Training Required: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No If no, explain:
If no, explain: No corrective action is required for this event as there are no changes to the current process, consumption of material, or how the product is produced. This added inspection guidelines are to avoid incorrect extrusion assembly defects.		

Title	Approval Name	Approval Signature	Date
Mgr. Quality Engineering	Hai Nguyen		25 JUL 2023
Mgr. Manufacturing Engineering	Jake Stanislowski		25 JUL 2023
Mgr. Operations	Matthew Benson		25 JUL 2023

CONTROLLED COPY

Deviation Authorization

FM0002.RevF

① UK55, 23JW 2023



DA | 2484
2468
①

Description/Objectives of Training:
DA- Inspection at final QC, Op#1050.

Group Training Record

Procedure:

- 100% inspection at Op#1050 per the instructions below.
- Inspect 1 part at a time.
- Inspection is focused on the correct MM0179-01 and MM1536-01 assembly.
- Use the example MM0179-01 and MM1536-02 fixture for inspection. (See image 1)
① **MM0179-01** **MM1536-01** **type correction TS12** **10AUG23**

CONTROLLED COPY




Image- 1

Step 1:

- Visually locate the MM0180-01 (Vestamid) transition to MM0179-01 on the completed part approximately 9.75" from the distal end using magnification light 2.25X minimum.
- Align the fixture MM0179-01 extrusion proximal end to the Vestamid transition on completed part. (See image 2)



Image- 2

- Visually verify the MM0179-01 distal end of the fixture is approximately at the same location on the completed part. (See image 3)





Image- 3

- Scrap the part if the transition is not approximately aligned. Save the scrapped parts for Engineer review.
- If the part transition is aligned, move to Step 2.

Step 2:

- Visually verify the MM1536-01 distal end of the fixture is approximately at the same location on the completed part. (See image 4)



- Image- 4
- Scrap the part if the transition is not approximately aligned. Save the scrapped parts for Engineer review.
 - If the part transition is aligned, the part passes inspection.
 - Use Image 5 as a guide for GOOD and BAD extrusion transition alignment.

1	MM0179-01 GOOD PART	MM1536-01
2	MM1536-01	MM0179-01 MM0179-01 and MM1536-01 Wrong Order - BAD PART
3	MM0179-01	MM0179-01 Two MM0179-01 - BAD PART
4	MM1536-01	MM1536-01 Two MM1536-01 - BAD PART

Image - 5

Edited to 11/15/2023
Edited to 13 Feb 2024 5228 V9beta4

CONTROLLED COPY DEVIATION AUTHORIZATION NUMBER: DA2564

CREGANNA
MEDICAL
is part of



DEVIATION AUTHORIZATION FORM

Requestor Name: Krishna Selvaraj			
Document Number Affected	Revision		
Doc #3005206 (MPI0238)	BP		
Deviation From:	Deviation To:		
Doc #3005206 (Flex Commander MPI0238): OPER850.11: Using a laser micrometer, check the D1M06 outer diameter. Position the laser indicator as close to the distal edge as possible. Start the measurement, then slowly move the part through the laser micrometer until reaching the lower edge of the shoulder.	Doc #3005206 (Flex Commander MPI0238): OPER850.11: Using a laser micrometer at OPER900 (TMI0700-01) , check the D1M06 outer diameter. Position the laser indicator as close to the distal edge as possible. Start the measurement, then slowly move the part through the laser micrometer until reaching the lower edge of the shoulder.		
Justification: TMI0602 lasermic which is currently used in SA0155-01 Flex commander product at OPER850 for Dim 6 inspection has mechanical failure and confirmed as not usable. TMI0700-01 lasermic is used at OPER900 for 100% inspection for Dim 1, Dim 6 and Dim 9. Since TMI0700-01 is already qualified to inspect Dim 6 per ES0647: Laser micrometer equivalency test, there is no additional risk in using TMI0700-01 for OPER850 Dim 6 inspection till TMI0602 issue is resolved.			
Part Number Affected	Revision		
SA0155-01	H		
Start Date:	End Date:	Lot Number:	
16 Nov 23	15 DEC 23	N/A	
Risk Assessment: Is there any potential risk(s) that may occur as a result of the proposed deviation including the following: Control Plans <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No FMEA's <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Validations <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Details (if any): N/A			
If yes to any of the above, what controls are being put in place to mitigate the risk – N/A			
Corrective Action Required:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No		
If no, explain: This is a temporary change to use TMI0700-01. DA will be removed once the lasermic TMI0602 issues are resolved and accepted for usage.			
Training Required:	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	If no, explain: N/A	
Title	Approval Name	Approval Signature	Date
Engineering Manager	Jake Stanislowski		16 Nov 2023
Quality Manager	Jay Zabel		16 Nov 2023
Operations Manager	Matthew Benson		16 Nov 2023



Document No: 5105589
FM5104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 500000297999

OP 400

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
TM10745	44	12:30 Am	430	SY47	26 Jan 24	12:42 Am	415	SY47	26 Jan 24	16
TM10745	44	1:07 Am	428	SH85	26 Jan 24	1:19 Am	415	SH85	26 Jan 24	16
TM10745	44	1:20 Am	429	V078	26 Jan 24	1:32 Am	415	V078	26 Jan 24	16
TM10745	44	6:50 AM	430	AM47	26 JAN 24	7:02 AM	415	AM47	26 JAN 24	16
TM10745	44	7:20 am	430	OS21	26 Jan 24	7:32 am	415	OS21	26 Jan 24	16
TM10745	44	7:49 Am	430	SN67	26 Jan 24	8:01 Am	415	SN67	26 Jan 24	16
TM10745	44	8:20am	430	OS21	26 Jan 24	8:32am	415	OS21	26 Jan 24	16
TM10745	44	9:35am	430	OS21	26 Jan 24	9:47am	415	OS21	26 Jan 24	16
TM10745	44	10:10am	430	OS21	26 Jan 24	10:22am	415	OS21	26 Jan 24	16
TM10745	44	10:50am	430	OS21	26 Jan 24	11:02am	415	OS21	26 Jan 24	16
TM10745	44	11:22am	430	RV16	26 Jan 24	11:34am	415	RV16	26 Jan 24	16
TM10745	44	11:55am	430	AL67	26 Jan 24	12:07pm	415	AL67	26 Jan 24	16

① CB58 26Jan24

② 6155 29Jan24



Document No: 5105589
FM5104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 50000029799

OP 400



Document No: 5105589
FM5104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 500000297999

OP 400

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
Tm10942	44	12:50 Am	430	SY47	26Jan24	1:02 Am	415	SY47	26Jan24	16
Tm10942	44	1:20 Am	429	SH85	26Jan24	1:32 Am	415	SH85	26Jan24	16
Tm10942	44	1:55 Am	430	SA85	26Jan24	2:07 Am	415	V078	26Jan24	16
Tm10942	44	6:35 Am	430	OS21	26Jan24	6:47 am	415	OS21	26Jan24	16
Tm10942	44	7:05 am	430	OS21	26Jan24	7:17 am	415	OS21	26Jan24	16
Tm10942	44	7:35 Am	430	SN67	26Jan24	7:47 Am	415	SN67	26Jan24	16
Tm10942	44	8:05am	429	OS21	26Jan24	8:17am	415	OS21	26Jan24	16
Tm10942	44	8:35am	430	OS21	26Jan24	8:47am	415	OS21	26Jan24	16
Tm10942	44	9:55am	430	OS21	26Jan24	10:17am	415	OS21	26Jan24	16
Tm10942	44	10:25am	430	OS21	26Jan24	10:37am	415	OS21	26Jan24	16
Tm10942	44	11:10am	430	OS21	26Jan24	11:22am	415	OS21	26Jan24	16
Tm10942	44	11:40Am	430	CD19	26Jan24	11:52Am	415	CD19	26Jan24	16



Document No: 5105589

FM5104665 Rev: C

Document Type: Manufacturing Form

Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 50000297999

OP 400



Document No: 5106073

Rev: E

Document Type: Manufacturing Form

Title: SA0155-01 Visual Rework Form

PO #: 500000297999OP #: 500 Shift #: 3rd**Total Parts Reworked:**162

Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	<u> </u>	<u>2</u>
EH	Exposed Hypotube	<u> + </u>	<u>45</u>
EW	Exposed Wire	<u> + </u>	<u>90</u>
MP	Micropores	<u>N/A</u>	<u>N/A</u>
SCR	Scratch	<u> </u>	<u>5</u>
SKV	Skive Marks	<u> </u>	<u>8</u>
VD	Voids	<u> </u>	<u>12</u>
N/A	N/A	<u>N/A</u>	<u>N/A</u>
Inspected By (Sign and Date):		<u>LL61.2546.</u>	<u>26 Jan 24</u>

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only.Data Uploaded for Engineering Review (Check): ☒



Document No: 5106073
Rev: E
Document Type: Manufacturing Form
Title: SA0155-01 Visual Rework Form

PO #: 50000291999

OP #: 500 Shift #: 2nd

Total Parts Reworked:		<u>3</u>	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	<u>N/A</u>	<u>0</u>
EH	Exposed Hypotube	<u>/</u>	<u>1</u>
EW	Exposed Wire	<u> </u>	<u>2</u>
MP	Micropores	<u>N/A</u>	<u>0</u>
SCR	Scratch	<u>N/A</u>	<u>0</u>
SKV	Skive Marks	<u>N/A</u>	<u>0</u>
VD	Voids	<u>N/A</u>	<u>0</u>
<u>N/A</u>	<u>N/A</u>	<u>N/A</u>	<u>0</u>
Inspected By (Sign and Date):		<u>Anita</u> 25 Jan 24	

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only.

Data Uploaded for Engineering Review (Check):



Document No: 6102646

Rev: A

Document Type: Manufacturing Form

Title: SA0155-01 Tipping Rework Form

PO #: 506000297999

OP #: 750 Shift #: 3st

Total Parts Reworked:		90	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
DIM07 OS / WO	DIM07 Oversized (Window Open)		20
DIM07 US / WC	DIM07 Undersized (Window Closed)	~1A	~1A
EH	Exposed Hypotube	~1A	50
N/A	Glue stopper		20
Inspected By (Sign and Date):		DY29	26 Jan 24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):

PRODUCTION ORDER# 500000297999

OP 800

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
TM10409	N/A	7:14AM	190°F	PL22	26 Jan 24	8:24AM	190°F	PL22	26 Jan 24	27
TM102036	N/A	7:55AM	190°F	PL22	26 Jan 24	9:05AM	190°F	PL22	26 Jan 24	45
TM10409	N/A	8:34AM	190°F	PL22	26 Jan 24	9:44AM	190°F	PL22	26 Jan 24	42
TM102036	N/A	9:32AM	190°F	YKHO	26 Jan 24	10:42AM	190°F	YKHO	26 Jan 24	25
TM10409	N/A	10:08AM	190°F	YKHO	26 Jan 24	11:18AM	190°F	YKHO	26 Jan 24	41
TM102036	N/A	10:44AM	190°F	YKHO	26 Jan 24	11:54AM	190°F	YKHO	26 Jan 24	28
TM10409	N/A	11:56AM	190°F	YKHO	26 Jan 24	1:06PM	190°F	YKHO	26 Jan 24	55
TM10409	N/A	1:10PM	190°F	BD64	26 Jan 24	2:20PM	190°F	BD64	26 Jan 24	31
TM102036	N/A	1:50PM	190°F	BD64	26 Jan 24	3:10PM	190°F	BD64	26 Jan 24	30
TM10409	N/A	3:05PM	190°F	BD64	26 Jan 24	4:15PM	190°F	BD64	26 Jan 24	71
TM10409	N/A	4:25PM	190°F	BD64	26 Jan 24	5:35PM	190°F	BD64	26 Jan 24	27
TM10409	N/A	4:50AM	190°F	KL45	27 JAN 24	6:00 AM	190°F	KL45	27 JAN 24	63
N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A



Document No: 6102619

Rev: B

Document Type: Manufacturing Form

Title: SA0155-01 Dimensional/Visual Rework Form

PO #: 500000297999

OP #: 900 Shift #: 3

Total Parts Reworked:		135	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	XX XX //	12
EH	Exposed Hypotube	XX // //	9
EW	Exposed Wire	XX //	61
MP	Micropores	XX XX //	13
SCR	Scratch	XX //	82
SKV	Skive Marks	XX //	6
VD	Voids	XX XX XX XX XX //	26
DIM01 US	DIM01 OD Undersized	N/A	N/A
DIM06 US	DIM06 OD Undersized	XX XX //	13
DIM06 OS	DIM06 OD Oversized	XX //	6
DIM09 US	DIM09 OD Undersized	N/A	N/A
Inspected By (Sign and Date):		DL07 MC17 KX54	26 Jan 24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):

Maximum Force Reached During Tensile Test (10 samples accepted from final inspection for each lot shall be selected and tensile tested)																
Sample # -->	1	2	3	4	5	6	7	8	9	10	Avg	St Dev	K	Calculated Lower bound	Min Spec	Pass / Fail
Seg A	30.43	25.37	26.93	36.52	27.09	32.29	28.79	22.23	30.15	27.83	28.763	3.9233774	4.378	11.5864539	8.542	PASS
Seg B	78.96	80.27	76.73	69.89	70.44	72.12	78.95	74.3	77.39	66.34	74.539	4.6705043	3.981	55.9457226	8.542	PASS
Seg C	73.09	77.29	77.11	76.1	86.13	83.1	86.59	83.64	76.5	85.64	80.519	4.9875945	2.911	66.0001124	8.542	PASS

All Force Values are recorded in Pound-Force and Distance is in Inches
Specification for lower bound is 38N was converted to 8.542Lbf
First Peak Force was collected during test and has been included in the raw data file (this information will not be captured / summarized in the DA due to it is not required to used for DA acceptance.

EDW Commander Flex - Bend and Tensile Strength Testing


LOT #: 500000297999

Date: 29JAN2024

Inspector Name: AUGUSTINE JAH

Equipment ID: TMI0311B

Cal Due Date: 27 OCT 24



29 JAN 2024