

Production Order: 500000304680



Production Order Document
Production Order Qty: 500
PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Material Type: ZFRT Description: Edwards Flex Shaft Commander 155885
Production Version: 7987
Plant / Business Unit: 1213 / AC5

Order Type: ZSTD

Project Phase:

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
50	KITTING3 Kitting Devices 	Kitting Devices Perform Order Kitting, Load Minor Mandrels, Dry Extrusions, and Cut FEP Record Time Extrusions Enter Dryer (Initial/Time/Date): <u>GS85 9:30AM</u> 06 Feb 24 Record Time Extrusions First Exit Dryer (Initial/Time/Date): <u>K82 2:15pm</u> 07 Feb 24 Record Dryer Shelf #: <u>N/A</u>	N/A	N/A	05 FEB 24	BV57

Notes: DA 2484, 2564

N/A
N/A

Date Printed: 02/05/2024 / 19:20:32

Page: 1 of 18



SA0155-01

CREGANNA MEDICAL
is part of



Production Order: 500000304680



Production Order Document
Production Order Qty: 500
PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details					Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	RM0158-01	E	F	PC	200	N/A <u>88018</u>	N/A 200		
		1000-1153-01	A	A	PC	594	N/A <u>89090</u> <u>89092</u> <u>88694</u>	N/A 200 200 200		
		1000-2053-01	A	A	PC	500	<u>0000287543</u>	500		
		MM1537-02	A	A	PC	500	<u>0000288401</u>	500	N/A	N/A
		TL0167-02	E	E	PC	70	N/A	Bulk		
		TL0165-05	J	J	PC	5	N/A	Bulk		
		TL0165-03	J	J	PC	5	N/A	Bulk		
							N/A	Bulk		

Notes:

N/A

N/A

N/A

Date Printed: 02/05/2024 / 19:20:32

Page: 2 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000304680



Production Order Document
Production Order Qty: 500
PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details						Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	141967-01	02	02	PC	500	① 0839985794 85502	204				
	RM7349-02	C	C	PC	543	88396 82864	300				
	RM7348-01	C	C	PC	500	85673	600	N/A			
	RM4001-01	B	B	PC	125	88702	500	N/A			
	RM0607-01	D	D	PC	56	82817 78322	100	N/A	N/A	N/A	N/A
	RM0498-01	C	C	PC	500	0000287647	25	482			
	RM0009-04	I	I	PC	1	88992	N/A	N/A			
	RM0009-04	I	I	PC	1	88992	Bulk	Bulk			

Notes:

N/A

N/A

N/A

Date Printed: 02/05/2024 / 19:20:32 ① KPD2 08Feb24

Page: 3 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000304680



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details					Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	MM1538-01	A	A	PC	500	<u>N/A</u>	Bulk		
							<u>0000290562</u>	<u>500</u>		
		MM1537-01	A	A	PC	1000	<u>N/A</u>	<u>N/A</u>		
							<u>0000290561</u>	<u>1120</u>		
		MM0177-01	C	C	PC	500	<u>N/A</u>	<u>N/A</u>		
							<u>0000284208</u>	<u>500</u>		
		MM0180-01	E	E	PC	500	<u>N/A</u>	<u>N/A</u>		
							<u>0000287541</u>	<u>500</u>		
		MM0178-01	E	E	PC	500	<u>N/A</u>	<u>N/A</u>		
							<u>0000276174</u>	<u>500</u>		
		MM0176-01	D	D	PC	500	<u>00002711050</u>	<u>80</u>		
							<u>0000288413</u>	<u>500</u>		
		MM0074-01	G	G	PC	500	<u>N/A</u>	<u>N/A</u>		
							<u>0000300897</u>	<u>521</u>		
							<u>0000292833</u>			
							<u>① 0000228</u>	<u>30</u>		

Notes:

N/A

N/A

N/A

Date Printed: 02/05/2024 / 19:20:32 ① KPO2 08 Feb 24

Page: 4 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000304680Production Order Document
Production Order Qty: 500**Material: SA0155-01 Rev F**PC
Sheet: 1 of 1

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	N/A	N/A	N/A	N/A	N/A
100	CATASY01 Catheter Assembly 1 	Line Clearance Perform Line Clearance and Heat Gun Setting	500	0	07Feb24	V078
	Line Clearance					
	Confirmation Reqd(Milestone)					
150	CATASY01 Catheter Assembly 1 	Major and Minor Mandrel Assembly	500	0	07Feb24	JY90 CL30 TRN SH85 SH23
	Major and Minor Mandrel Assembly					
Notes:						
N/A						
N/A						
N/A						

Date Printed: 02/05/2024 / 19:20:32

Page: 5 of 18



SA0155-01

CREGANNA
MEDICAL
is part of

Production Order: 500000304680



Production Order Document
Production Order Qty: 500

PC

Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
NIA	Confirmation Reqd(Milestone)	NIA	NIA	NIA	NIA	NIA
200	CATASY01 Catheter Assembly 1 	Loading Braid Stock	500	0	07Feb24	9531 7835 V078
	Loading Braid Stock					
	Confirmation Reqd(Milestone)					
250	CATASY01 Catheter Assembly 1 	Trim Braid Wire at Proximal End	500	0	07 Feb 24	C130
Notes:						
NIA						
NIA						
NIA						

Date Printed: 02/05/2024 / 19:20:32

Page: 6 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000304680



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	Trim Braid Wire at Proximal End N/A Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
300	CATASY01 Catheter Assembly 1 Insert Cut Hypo Tube Insert Cut Hypo Tube Confirmation Reqd(Milestone)	Insert Cut Hypo Tube	500	0	07FEB24	ST93 ST96
350	CATASY01 Catheter Assembly 1	Load Tubing	500	0	07FEB24	CLO5 DX35 GS22
Notes: N/A N/A N/A						

Date Printed: 02/05/2024 / 19:20:32

Page: 7 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000304680



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	NIA Load Tubing Confirmation Reqd(Milestone)	NIA	NIA	NIA	NIA	NIA
400	CATASY01 Catheter Assembly 1 Reflow Confirmation Reqd(Milestone)	Reflow	500	0	07Feb24 CL30 N078 S#85	
450	CATASY01 Catheter	FEP Removal	500	0	07Feb24 S790 Y014	
Notes: NIA NIA NIA						

Date Printed: 02/05/2024 / 19:20:32

Page: 8 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000304680



Production Order Document
Production Order Qty: 500

Material: SA0155-01 Rev F

PC
Sheet: 1 of 1

Cpr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	Assembly 1 FEP Removal N/A Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
500	CATASY01 Catheter Assembly 1 In-process Inspection and Rework Confirmation Reqd(Milestone)	In-process Inspection and Rework Material Consumed: Part #: 1000-1153-0 Batch #: 88747 Qty: N/A Part #: N/A Batch #: N/A Qty: N/A	482	EW - HHT HHT " DF - IIII VD - I SKV - I 18	P266 VL91 PP40 FD415 07Feb24	
N/A	N/A	N/A N/A N/A N/A	N/A	N/A	N/A	N/A
Notes:						

Date Printed: 02/05/2024 / 19:20:32

Page: 9 of 18



SA0155-01

CREGANNA
MEDICAL
is part of





Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
550	CATASY01 Catheter Assembly 1 	Remove Heat Shrink & Mandrel Remove Heat Shrink & Mandrel Confirmation Reqd(Milestone)	482	0	07Feb24	Y936 MV78 RS23 FBO1
600	CATASY01 Catheter Assembly 1 	Distal Tip Assembly Distal Tip Assembly Confirmation	475	DL-IH1 MAS-1 ⑤① 08Feb24 ⑦	MM02 MR60 VA96	

Notes:

MA

MA

MA

Date Printed: 02/05/2024 / 19:20:32

Page: 10 of 18



SA0155-01

CREGANNA MEDICAL
is part of

08446 08 Feb 24



Material: SA0155-01 Rev F

PC
Sheet: 1 of 1

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
NP	Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
650	CATASY01 Catheter Assembly 1 Loading Heat Shrink Confirmation Reqd(Milestone)	Loading Heat Shrink	475	0 08Feb24	ML38 D429	
700	CATASY01 Catheter Assembly 1 Tipping	Tipping Record Tipping Oven Information: TMI: 0521 Cal Due: 31 May 24 TMI: 2083C Cal Due: 31 May 24 TMI: 0356 Cal Due: 31 May 24 TMI: 0936A Cal Due: 31 May 24	475	0 08Feb24	ML38 RS23	

Notes:

N/A
N/A
N/A

Date Printed: 02/05/2024 / 19:20:32



SA0155-01



Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
750	CATASY01 Catheter Assembly 1 	<p>Tip Inspection/ Flash Removal Material Consumed:</p> <p>Part #: Pm0155-01 Batch #: 88702 Qty: 10 Part #: pm060201 Batch #: 78322 Qty: 4 Part #: N/A Batch #: N/A Qty: N/A Part #: N/A Batch #: N/A Qty: N/A Part #: N/A Batch #: N/A Qty: N/A</p>	475	0	08 Feb 24	SV46 PH59 HV36
800	CATASY01 Catheter Assembly 1 	Major Mandrel Removal	H66	AED-#111 ①	④ 08 Feb 24	SG88 XL91 SSH4

Notes:

Date Printed: 02/05/2024 / 19:20:32

Page: 12 of 18



SA0155-01

④ SSH4 08 Feb 24

CREGANNA
MEDICAL
is part of



Material: SA0155-01 Rev F

PC
Sheet: 1 of 1

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	Major Mandrel Removal N/A Confirmation Reqd(Milestone)		N/A	N/A	N/A	N/A
850	CATASY01 Catheter Assembly 1 	Cut to Length Record DIM05 gage result for the first 5 parts at the start of operation: 1. <u>passed</u> 2. <u>passed</u> 3. <u>passed</u> 4. <u>passed</u> 5. <u>passed</u>	464	SKV-1 EW-1 (2)	08 Feb 24	MILS5552
900	QUALITY1 Quality Inspection & Review	Quality Inspection and Review Perform Quality Inspection per QIP Document #3107610 Record Data in SAP ROS	N/A	N/A	N/A	HT72 MV33 SH04

Notes:

	N/A
	N/A
	N/A

Date Printed: 02/05/2024 / 19:20:32



SA0155-01



Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
	 Quality Inspection & Review Confirmation Reqd(Milestone) <i>N/A</i>	<p>Re-Inspect after re-work.</p> <p>Required Inspection Visual/OD Inspection Record Inspection Data in SAP ROS Record Laser Micrometer Information:</p> <p>TMI: <u>0700-01</u> Cal Due: <u>31 May 24</u> TMI: <u>N/A</u> Cal Due: <u>N/A</u> Material Consumed: Part #: <u>RM4001-01</u> Batch #: <u>88702</u> Qty: <u>N/A</u> Part #: <u>RM0607-01</u> Batch #: <u>78322</u> Qty: <u>N/A</u> Part #: <u>1000-1153-01</u> Batch #: <u>88694</u> Qty: <u>N/A</u> Part #: <u>RM0158-01</u> Batch #: <u>88018</u> Qty: <u>N/A</u> Part #: <u>N/A</u> Batch #: <u>N/A</u> Qty: <u>N/A</u> </p>	419	MAR-HHHH HT DIS-HTHTHT III DEL-II #605-HTI #705-I DL-III #805-I #805-I EW-II #905-I (45)	XL91 KL67 K155 P146 KT47 <i>Open 2/24</i>	
950	QUALITY1 Quality Inspection & Review	<p>Quality Inspection & Review Borescope Inspection Record Inspection Data in SAP ROS Record Tip Gage Information: TMI: <u>50713B</u> Cal Due: <u>12 Apr 24</u> Record Caliper Information:</p>	N/A	N/A	N/A	N/A

Notes:

Date Printed: 02/05/2024 / 19:20:32

Page: 14 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Material: SA0155-01 Rev F

PC
Sheet: 1 of 1

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	Quality Inspection & Review Confirmation Reqd(Milestone)	TMI: 0733 Cal Due: 30 Apr 24 Record DIM02 Go/No-Go Gage Information: TMI: 0691 Cal Due: 30 Sep 25 TMI: 0692 Cal Due: 30 Sep 25 Record DIM02 Inspection Results N = 54: Pass: 54 Fail: 0	386	DIS-HH HT HT HT HT STR-11 ① DIS-HH SD WK-1 DIS-LL (33)	08 Feb 24	XL91 KL67 0521
1000	QUALITY1 Quality Inspection & Review Quality Inspection & Review Confirmation Reqd(Milestone)	Quality Inspection & Review Leak Test Record Inspection Data in SAP ROS Record Leak Tester Information: TMI: 1056 Cal Due: 31 May 24 Record Length Gage Information: TMI: 0889D Cal Due: 30 Sep 24 Record Calibrated Ruler Information: TMI: 0629 Cal Due: 30 Sep 24	370	LT-HH HT HT 1 (16)	08 Feb 24	XL91 KL67 SS44

Notes:
N/A
N/A
N/A


Date Printed: 02/05/2024 / 19:20:32



SA0155-01

(OP146 08 Feb 24)

CREGANNA
MEDICAL
is part of





Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp Qty.	Scrap Qty & Desc.	Date Comp.	Initials
N/A	N/A	N/A	N/A	N/A	N/A	N/A
1050	QUALITY1 Quality Inspection & Review  Quality Inspection & Review Confirmation Reqd(Milestone) 	Required Inspection Visual Final Inspection Perform Quality Inspection per QIP Document #3107610 Record Data in SAP ROS	349	SCR - HII (TT) Del - III (TT) FM - II (TT) SCR - H SCR - I DIS - IIII VD - IIII Mar - II	(21)	SV43 08Feb24
1100	CATASY01 Catheter Assembly 1  Line Closure	Line Closure Perform Line Closure Settle materials issued to production order (Initials/Date): <u>KP02 08Feb24</u>	N/A	N/A	08Feb24	KP02

Notes:

 N/A
 N/A
 N/A

Date Printed: 02/05/2024 / 19:20:32

① SV43 08 Feb 24 Page: 16 of 18



SA0155-01

 CREGANNA
 MEDICAL
 is part of


Production Order: 500000304680



Production Order Document
Production Order Qty: 500
PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Opr No.	Planned WorkCenter Description	Operation Details	Comp. Qty.	Scrap Qty & Desc.	Date Comp.	Initials
1100	Confirmation Reqd(Milestone)	N/A	N/A	N/A	N/A	N/A
1150	PACKINT1 Packing assembly Package Confirmation Reqd(Milestone)	Package Package, Label, and Ship Finished Parts	349	0 09 Feb 24	AP 10 AP 10 09 Feb 24	

Notes:

N/A AP 10 09 Feb 24

Date Printed: 02/05/2024 / 19:20:32

Page: 17 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Production Order: 500000304680



Production Order Document
Production Order Qty: 500

PC
Sheet: 1 of 1

Material: SA0155-01 Rev F

Batch Number: 0000304680

By: AP10

Date: 09 Feb 24

Reviewed By:

RB29

Date:

09 feb 24

Notes:

N/A AP10 09 Feb 24

Date Printed: 02/05/2024 / 19:20:32

Page: 18 of 18



SA0155-01

CREGANNA
MEDICAL
is part of



Basis to 2024 3228 11/16/23
Excess to 19 Feb 2024 3228 11/16/23
CREGANNA MEDICAL is part of



CREGANNA MEDICAL
is part of

Extend to 2023 3228 11/16/23
Basis to 2023 3228 11/16/23


Requestor Name: Udhesh Kapadnis

DEVIATION AUTHORIZATION NUMBER: 2484
* See attached email extension to 2484
TS12 24 AUG 23 3228 11/16/23

CONTROLLED COPY

Deviation From:	Deviation To:	
QIP3107610, Section 8.0 Inspection Requirements (Supplemental Visual Inspection) OP 1050: Current QIP3107610 does not state to inspect for the correct extrusion configuration.	<p>This DA allows addition inspection for correct assembly of extrusion material MM0179-01 and MM1536-01 during performing QIP3107610, Section 8.0 Inspection Requirements (Supplemental Visual Inspection) OP 1050.</p> <p>See instructions attached to this DA.</p> <p>Pre-Release Assembly Diagram</p>	
Document Number Affected	Revision	
3107610	L	

Risk Assessment: Is there any potential risk(s) that may occur as a result of the proposed deviation including the following: Control Plans <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No FMEAs <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Validations <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Details (if any): N/A	Part Number Affected <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">SA0155-01</td> <td style="width: 50%;">H</td> </tr> </table> Start Date: 26 Jul 2023 End Date: 25 Aug 2023 Lot Number: N/A			SA0155-01	H
SA0155-01	H				
Corrective Action Required: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No If no, explain: _____					
If no, explain: No corrective action is required for this event as there are no changes to the current process, consumption of material, or how the product is produced. This added inspection guidelines are to avoid incorrect extrusion assembly defects.					
Training Required: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No If no, explain: _____					

Title	Approval Name	Approval Signature	Date
Mgr. Quality Engineering	Hai Nguyen		25 JUL 2023
Mgr. Manufacturing Engineering	Jake Stanislawski		25 JUL 2023
Mgr. Operations	Matthew Benson		25 JUL 2023

FM0002.RevF Deviation Authorization

CONTROLLED COPY

① UK55, 23JW 2023



DA | DA | 2484
2468 • ①

Description/Objectives of Training:
DA- Inspection at final QC, Op#1050.

Group Training Record

Procedure:

- 100% inspection at Op#1050 per the instructions below.
- Inspect 1 part at a time.
- Inspection is focused on the correct MM0179-01 and MM1536-01 assembly.
- Use the example MM0179-01 and MM1536-02 fixture for inspection. (See image 1)
MM0179-01 type connection TS12 10AUG23

CONTROLLED COPY




Image- 1

Step 1:

- Visually locate the MM0180-01 (Vestamid) transition to MM0179-01 on the completed part approximately 9.75" from the distal end using magnification light 2.25X minimum.
- Align the fixture MM0179-01 extrusion proximal end to the Vestamid transition on completed part. (See image 2)



Image- 2

- Visually verify the MM0179-01 distal end of the fixture is approximately at the same location on the completed part. (See image 3)




Image- 3

- Scrap the part if the transition is not approximately aligned. Save the scrapped parts for Engineer review.
- If the part transition is aligned, move to Step 2.

Step 2:

- Visually verify the MM1536-01 distal end of the fixture is approximately at the same location on the completed part. (See image 4)




Image- 4

- Scrap the part if the transition is not approximately aligned. Save the scrapped parts for Engineer review.
- If the part transition is aligned, the part passes inspection.
- Use Image 5 as a guide for GOOD and BAD extrusion transition alignment.

1	MM0179-01 GOOD PART	MM1536-01
2	MM1536-01	MM0179-01 MM0179-01 and MM1536-01 Wrong Order - BAD PART
3	MM0179-01	MM0179-01 Two MM0179-01 - BAD PART
4	MM1536-01	MM1536-01 Two MM1536-01 - BAD PART

Image - 5

Edo to HENRY J228 11/16/03

Edo to 13 Feb 2024 5228 1/1/04

CONTROLLED COPY DEVIATION AUTHORIZATION NUMBER: DA2564

CREGANNA
MEDICAL
is part of



DEVIATION AUTHORIZATION FORM

Requestor Name: Krishna Selvaraj			
Document Number Affected	Revision		
Doc #3005206 (MPI0238)	BP		
Deviation From:	Deviation To:		
Doc #3005206 (Flex Commander MPI0238): OPER850.11: Using a laser micrometer, check the DIM06 outer diameter. Position the laser indicator as close to the distal edge as possible. Start the measurement, then slowly move the part through the laser micrometer until reaching the lower edge of the shoulder.	Doc #3005206 (Flex Commander MPI0238): OPER850.11: Using a laser micrometer at OPER900 (TMI0700-01) , check the DIM06 outer diameter. Position the laser indicator as close to the distal edge as possible. Start the measurement, then slowly move the part through the laser micrometer until reaching the lower edge of the shoulder.		
Justification: TMI0602 lasermic which is currently used in SA0155-01 Flex commander product at OPER850 for Dim 6 inspection has mechanical failure and confirmed as not usable. TMI0700-01 lasermic is used at OPER900 for 100% inspection for Dim 1, Dim 6 and Dim 9. Since TMI0700-01 is already qualified to inspect Dim 6 per ES0647; Laser micrometer equivalency test, there is no additional risk in using TMI0700-01 for OPER850 Dim 6 inspection till TMI0602 issue is resolved.			
Part Number Affected	Revision		
SA0155-01	H		
Start Date:	End Date:		
16 Nov 23	15 DEC 23		
Lot Number:	N/A		
Risk Assessment: Is there any potential risk(s) that may occur as a result of the proposed deviation including the following: Control Plans <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No FMEA's <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Validations <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Details (if any): N/A			
If yes to any of the above, what controls are being put in place to mitigate the risk - N/A			
Corrective Action Required: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No			
If no, explain: This is a temporary change to use TMI0700-01. DA will be removed once the lasermic TMI0602 issues are resolved and accepted for usage.			
Training Required: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No If no, explain: N/A			
Title	Approval Name	Approval Signature	Date
Engineering Manager	Jake Stanislowski		16 Nov 2023
Quality Manager	Jay Zabel		16 Nov 2023
Operations Manager	Matthew Benson		16 Nov 2023



Document No: 5105589
FM15104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 500000304680

OP 400



Document No: 5105589
FM5104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER#: 50000304680

OP 400

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
Tm10942	44	2:48pm	430	0521	07Feb24	3:00pm	415	0521	07Feb24	16
Tm10942	44	4:25pm	430	JY90	07Feb24	4:37pm	415	JY90	07Feb24	16
Tm10942	44	5:02pm	430	Sy47	07Feb24	5:14pm	415	Sy47	07Feb24	16
Tm10942	44	5:25pm	427	SH85	07Feb24	5:37pm	415	SA07	07Feb24	16
Tm10942	44	6:30pm	430	SH85	07Feb24	6:42pm	415	SH85	07Feb24	16
Tm10942	44	6:59pm	430	SH85	07Feb24	7:11pm	415	SH85	07Feb24	16
Tm10942	44	7:25pm	428	SH85	07Feb24	7:37pm	415	SH85	07Feb24	16
Tm10942	44	7:56pm	429	SH85	07Feb24	8:08pm	415	V078	07Feb24	16
Tm10942	44	8:52pm	430	V078	07Feb24	9:04pm	415	Sy47	07Feb24	16
Tm10942	44	9:35pm	430	V078	07Feb24	9:47pm	415	Sy47	07Feb24	16
Tm10942	44	10:25pm	430	V078	07Feb24	10:37pm	415	V078	07Feb24	16
Tm10942	44	10:55pm	430	V078	07Feb24	11:07pm	415	V078	07Feb24	16



Document No: 5105589
FM104665 Rev: C
Document Type: Manufacturing Form
Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 500000304680

OP 400



Document No: 5105589

FM5104665 Rev: C

Document Type: Manufacturing Form

Title: SA0155-01 Reflow Log Sheet Form

PRODUCTION ORDER# 500000304680

OP 400

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
Tm10745	44	4:10 pm	430	JY90	07Feb24	4:22 pm	415	JY90	07Feb24	16
Tm10745	44	4:41 pm	430	SH85	07 Feb 24	4:53 pm	415	SH85	07 Feb 24	16
Tm10745	44	5:17 pm	429	Sy47	07Feb24	5:29 pm	415	Sy47	07Feb24	16
Tm10745	44	6:20 pm	430	CL30	07Feb24	6:32 pm	415	CL30	07Feb24	16
Tm10745	44	6:44 pm	428	SH85	07Feb24	6:56 pm	415	SH85	07Feb24	16
Tm10745	44	7:14 pm	430	CL30	07Feb24	7:26 pm	415	CL30	07Feb24	16
Tm10745	44	7:38 pm	428	SH85	07Feb24	7:50 pm	415	SH85	07Feb24	16
Tm10745	44	9:07 pm	430	V078	07 Feb 24	9:19 pm	415	Sy47	07Feb24	16
Tm10745	44	9:25 pm	429	Sy47	07Feb24	9:37 pm	415	Sy47	07Feb24	16
Tm10745	44	9:50 pm	⑥ 429	Sy47	07Feb24	10:02 pm	415	Sy47	07Feb24	16
Tm10745	44	10:14 pm	429	V078	07Feb24	10:26 pm	415	V078	07Feb24	16
Tm10745	44	10:40 pm	429	V078	07Feb24	10:52 pm	415	V078	07Feb24	16

① Sy47 07Feb24




Document No: 5106073
Rev: E
Document Type: Manufacturing Form
Title: SA0155-01 Visual Rework Form

PO #: 50000304680 OP #: 500 Shift #: 2

Total Parts Reworked:		29	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	N/A	0
EH	Exposed Hypotube		4
EW	Exposed Wire		25
MP	Micropores	N/A	0
SCR	Scratch	N/A	0
SKV	Skive Marks	N/A	0
VD	Voids		4
N/A	N/A	N/A	0
Inspected By (Sign and Date):		Craig	07 Feb 24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only.

Data Uploaded for Engineering Review (Check):





Document No: 5106073

Rev: E

Document Type: Manufacturing Form


Title: SA0155-01 Visual Rework Form

PO #: 50000304680 OP #: 500 Shift #: 2nd

Total Parts Reworked:		24	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	n/a	n/a
EH	Exposed Hypotube	n/a	n/a
EW	Exposed Wire		17
MP	Micropores	n/a	n/a
SCR	Scratch	///	3
SKV	Skive Marks	n/a	n/a
VD	Voids		7
n/a	n/a	n/a	n/a
Inspected By (Sign and Date):		Varunegi Lor 07 Feb 24	

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only.

Data Uploaded for Engineering Review (Check):





Document No: 5106073

Rev: E

Document Type: Manufacturing Form


Title: SA0155-01 Visual Rework Form

PO #: 300000304680OP #: 500 Shift #: 2nd

Total Parts Reworked:		25	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	N/A	N/A
EH	Exposed Hypotube		1
EW	Exposed Wire		18
MP	Micropores		
SCR	Scratch	N/A	
SKV	Skive Marks	PP40 07 feb 24	
VD	Voids		6
N/A	N/A	N/A	N/A
Inspected By (Sign and Date):		PP40 07 feb 24	

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only.

Data Uploaded for Engineering Review (Check):





Document No: 6102646

Rev: A

Document Type: Manufacturing Form

Title: SA0155-01 Tipping Rework Form


PO #: 500000304680

OP #: 750 Shift #: 2nd

Total Parts Reworked:		40	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
DIM07 OS / WO	DIM07 Oversized (Window Open)		3
DIM07 US / WC	DIM07 Undersized (Window Closed)		28
EH	Exposed Hypotube		9
N/A	N/A	N/A	N/A
Inspected By (Sign and Date):		SV46 07 Feb 24	

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):



PO #: 500000304680OP #: 750 Shift #: 1st

Document No: 6102646

Rev: A


Document Type: Manufacturing Form

Title: SA0155-01 Tipping Rework Form

Total Parts Reworked:		<u>15</u>	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
DIM07 OS / WO	DIM07 Oversized (Window Open)	<u>N/A</u>	<u>N/A</u>
DIM07 US / WC	DIM07 Undersized (Window Closed)	<u> </u>	<u>5</u>
EH	Exposed Hypotube	<u> </u>	<u>3</u>
<u>N/A</u>	<u>Glue , Stopper</u>	<u> </u>	<u>7</u>
Inspected By (Sign and Date):		<u>PH59</u>	<u>08 Feb 24</u>

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):



PRODUCTION ORDER# 500000304680

OP 800

Oven #	Cycle #	Time In	Temp. In (Actual)	Initials	Date	Time Out	Temp. Out (Actual)	Initials	Date	Qty
Tm12036	N/A	7:31pm	190°F	SG88	07 Feb 24	8:32PM	190°F	XL91	07 Feb 24	41
Tm10409	N/A	8:37pm	190°F	SG88	07 Feb 24	9:47pm	190°F	SG88	07 Feb 24	32
Tm12036	N/A	9:19pm	190°F	SG88	① 07 Feb 24	10:29pm	190°F	SG88	07 Feb 24	43
Tm10409	N/A	10:06pm	190°F	SG88	07 Feb 24	11:16pm	190°F	SG88	07 Feb 24	48
Tm12036	N/A	10:47pm	190°F	SG88	07 Feb 24	11:57pm	190°F	SG88	07 Feb 24	40
Tm10409	N/A	11:46pm	190°F	SG88	07 Feb 24	12:56am	190°F	SG88	08 Feb 24	44
Tm12036	N/A	12:24AM	190°F	SG88	08 Feb 24	1:34AM	190°F	SG88	08 Feb 24	37
Tm10409	N/A	12:57AM	190°F	SG88	08 Feb 24	2:10am	190°F	SG88	08 Feb 24	50
Tm10409	N/A	4:20am	190°F	K155	08 Feb 24	5:30am	190°F	K155	08 Feb 24	34
Tm10409	N/A	5:35am	190°F	SS44	08 Feb 24	6:45am	190°F	SS44	08 Feb 24	36
Tm12036	N/A	6:00am	190°F	SS44	08 Feb 24	7:10am	190°F	SS44	08 Feb 24	36
Tm10409	N/A	6:50am	190°F	SS44	08 Feb 24	8:00am	190°F	SS44	08 Feb 24	25
N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

① SG88 07 Feb 24

PO #: 50000304680OP #: 900 Shift #: 2nd

Document No: 6102619
Rev: B
Document Type: Manufacturing Form
Title: SA0155-01 Dimensional/Visual Rework Form

Total Parts Reworked:		38	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	N/A	0
EH	Exposed Hypotube		3
EW	Exposed Wire		5
MP	Micropores	N/A	0
SCR	Scratch		20
SKV	Skive Marks	N/A	0
VD	Voids		1
DIM01 US	DIM01 OD Undersized	N/A	0
DIM06 US	DIM06 OD Undersized		15
DIM06 OS	DIM06 OD Oversized		1
DIM09 US	DIM09 OD Undersized	N/A	0
Inspected By (Sign and Date):		See H 07 Feb 24	

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):



Document No: 6102619
Rev: B
Document Type: Manufacturing Form
Title: SA0155-01 Dimensional/Visual Rework Form

PO #: 500000304680 OP #: 900 Shift #: 2

Total Parts Reworked:		45	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	N/A	N/A
EH	Exposed Hypotube	11	2
EW	Exposed Wire		23
MP	Micropores	N/A	N/A
SCR	Scratch		42
SKV	Skive Marks	N/A	N/A
VD	Voids	11	2
DIM01 US	DIM01 OD Undersized		
DIM06 US	DIM06 OD Undersized		
DIM06 OS	DIM06 OD Oversized		
DIM09 US	DIM09 OD Undersized		
Inspected By (Sign and Date):		HT 72 07 Feb 24	

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):



PO #: 500000304680 OP #: 900 Shift #: 1ST

Document No: 6102619
Rev: B
Document Type: Manufacturing Form
Title: SA0155-01 Dimensional/Visual Rework Form

Total Parts Reworked:		46	
Router Code	Defect Failure Mode	Reworkable Defects (Tally)	Total Defects
AB	Air Bubbles	N/A	N/A
EH	Exposed Hypotube	11	2
EW	Exposed Wire	HHH HHH HHH 11	17
MP	Micropores	N/A	N/A
SCR	Scratch	HHH HHH HHH	15
SKV	Skive Marks	1	1
VD	Voids	HHH	5
DIM01 US	DIM01 OD Undersized		N/A N/A
DIM06 US	DIM06 OD Undersized	HHH 1	6
DIM06 OS	DIM06 OD Oversized		N/A N/A
DIM09 US	DIM09 OD Undersized		N/A N/A
Inspected By (Sign and Date):		biss KT47	08 Feb 24

Note: Indicate tally marks in groups of 5. Scrap is to be recorded on the SAP router; this form is for reworked parts only. DIM01 OS, DIM09 OS, Foreign Material, and Cracks are not reworkable per MPI0238.

Data Uploaded for Engineering Review (Check):

Maximum Force Reached During Tensile Test (10 samples accepted from final inspection for each lot shall be selected and tensile tested)																
Sample # -->	1	2	3	4	5	6	7	8	9	10	Avg	St Dev	K	Calculated Lower bound	Min Spec	Pass / Fail
Seg A	24.74	26.78	24.68	26.2	24	25.97	24.63	24.71	26.36	23.6	25.167	1.0786313	4.378	20.444752	8.542	PASS
Seg B	58.33	60.49	67.16	63.78	66.05	61.53	60.88	60.6	63.68	60.04	62.254	2.8166103	3.981	51.0410746	8.542	PASS
Seg C	81.7	71.79	77.9	80.31	77.75	80.43	81.18	71.96	76.4	74.41	77.383	3.6788103	2.911	66.6739831	8.542	PASS

All Force Values are recorded in Pound-Force and Distance is in Inches
Specification for lower bound is 38N was converted to 8.542Lbf
First Peak Force was collected during test and has been included in the raw data file (this information will not be captured / summarized in the DA due to it is not required to used for DA acceptance.

EDW Commander Flex - Bend and Tensile Strength Testing

LOT #: 500000304680

Date: 08 FEB 24

Inspector Name: LUKASU C. TSHISHIMBI

Equipment ID: TMIO311B

Cal Due Date: 27 OCT 24

08Feb24
